

Cutting Instruments

Experts in high-quality cutting instruments for all applications in orthopedic surgery, including spine, trauma, extremities, and large joints

Design | Develop | Manufacture | Package | Deliver

Why Choose **Avalign** Cutting Instruments?

- Close tolerance grinding of market-leading cutting geometries for drills, reamers, taps, and fixation pins
- Able to meet challenging design requirements with innovation and specialized products, processes, equipment, and tooling
- State-of-the-art grind, mill, and swiss machines with robotic feeds for high volume products
- Polish, finish, clean and passivate, laser etch ... all in one location
- Significant new capacity installed and available to support your next project



Your Trusted CDMO Partner

Avalign is one of the only CDMOs that can provide every aspect of an orthopedic implant system, including: Implants and Coatings | Fixation Devices | Complex Instruments | Cutting Instruments | German Specialty Instruments | Cases & Trays | Fully Kitted Implant Banks and Instrument Sets

98%

On Time
Delivery
Commitment

98%

Customer Lot
Acceptance
Commitment

25+

Years
Experience

8-18

Week Lead
Time

Contact Us



(855) 282-5446



www.avalign.com